

0. EXECUTIVE SUMMARY (WHY THIS PRODUCT EXISTS)

Families don't want raw health data.

They want **assurance**, **context**, and **early awareness** — without panic.

This product enables **safe, consent-first family visibility** into daily health and routine changes using wearable and behavioural signals, designed especially for **elderly parents living independently** and **adult children living elsewhere (often globally)**.

1. PRODUCT VISION & PHILOSOPHY

1.1 Vision

To build a **family health & routine awareness platform** that helps people understand **how their loved ones are doing day-to-day**, without:

- medicalising daily life
- creating anxiety or false alarms
- invading privacy
- exposing elderly users to coercion or misuse

The product is about **awareness, not diagnosis**.

1.2 Core Insight

Families fundamentally ask:

“Is today normal for them?”

Not:

- “What was their heart rate at 3:12 PM?”

1.3 Phase-1 Principles

1. Wearable-first for credibility
2. Behavioural signals for context, not surveillance
3. Relative comparisons, not medical thresholds
4. Deterministic, explainable logic (no ML)
5. Primary-controlled data sharing only
6. Consent always precedes payment
7. Paid-only model with free trial
8. Trust > growth > virality

2. USER ROLES & MENTAL MODEL

2.1 Canonical Roles

Primary Member

- Person whose data is tracked
- Typically elderly parent / dependent

- Owns all data
- Controls sharing & permissions
- May be low-tech

Secondary Member

- Family member who views another member's data
- Typically adult child (often outside India)
- May be the payer
- Cannot access anything without explicit consent

Every user is always a **Primary Member for their own data**

A user may also be a **Secondary Member for others**

2.2 Mental Model

"You always see yourself first.

You only see others if they explicitly share."

No:

- access requests
- search / discovery
- social graph

3. DATA SOURCES

3.1 Wearable Data (Trial + Paid Only)

- Source: Spike API
- Devices: Garmin, Fitbit, Samsung, Oura, etc.
- Cloud-based sync
- Cost incurred **per connected device**

Used for:

- Vitals
- Sleep
- Activity
- Training load

3.2 Smartphone Data (Supportive Only)

Used **only** for behaviour & context:

- Smartphone usage time (daily total)
- Outside time (derived)
- Last active location (context)

- Steps (fallback only)

4. METRIC TAXONOMY

4.1 Vitals (Wearable Only)

- Heart Rate
- HRV
- Stress
- SpO₂
- Respiration
- Sleep
- Body Temperature
- Blood Pressure
- ECG (if supported)
- Blood Glucose (if supported)
- VO₂ Max

4.2 Activity

- Steps
- Calories
- Workouts
- Training load

4.3 Behaviour (Derived)

- Routine Consistency
- Smartphone usage time
- Outside time
- Last active location (context only)

5. DERIVED LOGIC & BASELINES

5.1 Baseline Calculation

- Rolling baseline over **14–21 valid days**
- Day excluded if:
 - o wearable not worn
 - o sync gaps exceed threshold
 - o phone data missing
- Baseline recalculated daily

5.2 Comparison Language (Phase-1)

All metrics use **relative states only**:

- Higher than usual
- Same as usual
- Lower than usual

No colours.

No medical thresholds.

5.3 Outside Time (Derived)

Logic

- Identify home cluster
- Measure time outside cluster
- Compare with baseline

Output

- Relative state only

5.4 Routine Consistency Score

Inputs

- Smartphone usage
- Outside time
- Steps
- Sleep (if wearable connected)

Logic

- Normalize vs baseline
- Behaviour weighted higher than vitals
- Output stability band:
 - o Stable
 - o Slightly different
 - o Noticeably different

Explainable. Deterministic. No ML.

6. BUSINESS MODEL (PAID-ONLY)

6.1 Paid-Only + 7-Day Free Trial

There is **no perpetual free tier** in Phase-1.

Rationale

- Core value requires wearable

- Phone-only users misjudge product
- Infra + Terra costs are real
- Target users already pay for wearables

Freemium deferred to **future funded phase**.

6.2 Trial Details

- Duration: 7 days
- Full product access
- Payment method required upfront

Supported Payments

- Credit / Debit Cards (global)
- UPI (India)

UX Copy

“Try the full experience free for 7 days.
No charge during trial. Cancel anytime.”

6.3 Trial Expiry Behaviour

If continued

- Seamless conversion

If cancelled / lapsed

- Wearable sync paused
- Data read-only
- Copy:

“Health tracking is paused. Restart subscription to continue.”

No deletion. No surprise charges.

7. PRICING

- ₹249 / month:
 - o 1 wearable
 - o Full insights
 - o 2 family member sharing
- Monthly recurring

8. CONSENT & PAYMENT MODEL (FINAL)

Consent always comes before payment.

Payment only unlocks already-approved capability.

There is **no post-payment approval anywhere**.

9. UX FLOWS — PAYMENT & ONBOARDING

Flow A — Add Yourself (Self-Pay)

1. Profile → Add Device
2. “Connect your wearable — ₹249/month
Try free for 7 days”
3. Start trial
4. Connect device
5. Sync begins

Flow B — Add Family Member’s Device (Delegated Pay)

Primary

“To connect a wearable, a trusted family member needs to enable this.”

CTA: Ask family member to enable
(Implicit consent)

Secondary

“[Parent’s name] wants to connect a wearable.”

Start trial → subscription active

Primary connects device.

Flow C — Add Additional Family Member (only upto 2 family members)

1. Primary invites
2. Secondary accepts

Invite itself = consent.

10. CORE APP UX

Home (Members Overview)

- Self first
- Others only if shared
- Metric deviation counts
- Last updated time

Member Detail

- Daily NLP summary
- Last active location
- Metrics grouped & ordered by deviation

Metric Detail

- 24h / 7d / 30d

- Graph + explanation
- Non-medical guidance

11. PERMISSIONS MODEL (FINAL)

11.1 Default Permission (IMPORTANT)

☒ All access enabled by default

Rationale:

- Most families want full visibility
- Reduces setup friction
- Matches emotional intent of sharing

11.2 Permission Categories

Category	Access
All	Everything
Vitals	Health metrics
Activity	Steps, workouts
Behaviour	Routine, phone usage

Multi-select supported.

11.3 Permission Modification UX

Primary Screen

Manage access for *[Name]*

- ☒ All
- ☒ Vitals
- ☒ Activity
- ☒ Behaviour

- Toggle individual metrics allowed; Individual metrics listed below their specific category

Changes apply instantly.

12. UX EDGE CASES (IMPORTANT)

Case 1: Primary revokes access

- Secondary loses access instantly
- No payment refund
- Slot freed

Case 2: Subscription lapses

- Wearable sync pauses
- Sharing pauses
- Old Data remains visible (read-only)

Case 3: Device disconnected

- Show “No recent data”
- No inferred values

Case 4: Permission reduced

- Secondary sees:

“Some data is no longer shared.”

- No prompts
- No upsell

Case 5: Multiple Secondary Members

- Each evaluated independently
- Permissions per Secondary

13. GTM RATIONALE (EXPLICIT)**Why Paid-First Works**

- Wearable owners already pay
- Family health = high intent
- Paid users give better feedback
- Less noise, better PMF signal
- Controlled burn during bootstrap

Future Freemium

- Introduced post-funding
- Used as top-of-funnel only

14. NON-GOALS (PHASE-1)

- Medical diagnosis
- Emergency alerts
- ML predictions
- Continuous GPS
- Doctor marketplace
- Fear-based nudges

15. ONE-LINE PRODUCT SUMMARY

A wearable-first family health awareness platform that helps people understand daily health and routine changes of loved ones — safely, clearly, and without medicalisation.